



IntelliVue Anesthetic Gas Modules – G1 and G5

Anesthesia monitoring choice

Versatile IntelliVue Anesthetic Gas Modules measure the five most commonly used anesthetic gases, as well as N₂O and CO₂. All provide inspiration and expiration values and calculate MAC values for display on Philips IntelliVue patient monitors. The IntelliVue G1 gas module measures the single agent chosen by the clinician. Advanced O₂ measurements is optional with the G1 and included standard with the G5.

The IntelliVue G5 features automatic agent identification and mixed-agent measurement capability.

Real-time gas monitoring data displayed on the patient monitor


The IntelliVue gas modules provide clinically solid measurements:

- Numerics, waveforms, and trends:
 - Sevoflurane
 - Isoflurane
 - Halothane
 - Desflurane
 - Enflurane
 - Nitrous oxide (N₂O)
 - Carbon dioxide (CO₂)
 - Oxygen (O₂)*
- Airway respiration rate
- MAC values
- Agent identification and measurements

Additional features include:

- Quick warm-up time
- Rapid response time
- Individual alarm limits for each parameter
- Automatic zero calibration
- Lifetime calibration for easier maintenance
- Highly durable solid-state sensor design

Efficiently monitor all commonly administered gases and view measurements in context with other physiologic parameters.



* O₂ optional on the IntelliVue G1.


Please ask your sales representative for details on compatibility.

Versatile gas monitoring for a range of anesthesia settings

The compact footprint of IntelliVue gas modules makes them practical for use in lower-acuity settings. At the same time, advanced capabilities such as automatic agent ID and mixed-agent measurement are available for higher-acuity environments.

Durable, compact, and easy to maintain

The gas sensors in the IntelliVue G1 and G5 gas modules have no moving parts, so they are reliable and durable. Zero calibrations are completed automatically during warm-up. A wide range of mounting options makes it easy to build your optimal anesthesia workspace.



IntelliVue G1 and G5 anesthetic gas modules are compact, yet powerful.

	IntelliVue G1	IntelliVue G1 with O ₂	IntelliVue G5
Anesthetic agents	✓	✓	✓
N ₂ O and CO ₂	✓	✓	✓
O ₂	✗	✓	✓
Agent identification	Manual	Manual	Automatic
Agent measurement ID	Single	Single	Mixed
Supported IntelliVue patient monitors	MP20 – MP50	MP20 – MP50	MP40 and higher*

* MP20 and higher in the US and Canada only.

M1013A IntelliVue G1 Anesthetic Gas Module

M1019A IntelliVue G5 Anesthetic Gas Module

Philips Commitment to Measurement Technologies

Philips is committed to providing best-in-class standard clinical measurements as well as innovative measurements to support clinicians' decisions at the patient's side.

- Maintaining and advancing the performance of existing, widely used standard-of-care measurements
- Investing heavily in research, development, and clinical validation of new, innovative parameters and algorithms
- Working with strategic partners to integrate next-generation measurements and technologies
- Providing interfaces to more than 100 third-party specialty measurement devices through the Philips VueLink module



© 2006 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Medical Systems is part of Royal Philips Electronics

www.medical.philips.com
 medical@philips.com
 fax: +31 40 27 64 887

Printed in The Netherlands
 4522 962 10451/862 * APR 2006

Philips Medical Systems
 Global Information Center
 P.O. Box 1286
 5602 BG Eindhoven
 The Netherlands