



High performance mobile light

Phililux S50 mobile

PHILIPS



S50 mobile is a high performance mobile light very effective in all types of surgeries, minor or major.

A right choice where mobility and high performance is a priority.

S50 mobile - Truly a moving performance.

Dome	1 No.
Dome diameter	630mm (nominal)
Spring arm stroke	375mm \pm 50mm
Field size	200mm diameter (max)
Light output	75,000lux
Colour temperature	4200°K \pm 200°K
Halogen bulb	5 Nos, 12V, 50W each
Low voltage unit with CVT	1 No.
Input supply	230 V AC \pm 10%, 50Hz \pm 2%
Intensity control	in 4 steps
Castor wheel	4 Nos. (2 Nos. with brake)
Max height from floor	2200mm \pm 50mm
Trolley size	890mm x 600mm (nominal)

HO/Western regional office tel : +91 022 6691 2000, Fax : +91 022 6691 2474
 Eastern regional office tel : +91 033 24559645, Fax : +91 033 24559423
 Northern regional office tel : +91 011 42095036, Fax : +91 011 26959798
 Southern regional office tel : +91 044 66501000, Fax : +91 044 66501075

Equipments are subject to change without notice. All changes comply with regulations governing the manufacture of medical equipment.



Philips Electronics India Ltd.
 8568 5191 | *09/08



ISO 9001-2000 Certificate No:142856